

UnionMedico®  
We do best for your drug

Semi

# Auto-Injector

Intramuscular • Subcutaneous



Super Grip™  
Auto-Injector




Color of choice

# Revolutionary

Auto-Injector System  
90° Reusable Auto-Injector

100% Depot control  
Controlled delivery



## Easy - To - Use - Patients controlled injection


- Greater drug efficiency
- Placing Needle in tissue
- No more anxiety and needle phobia
- Right injection angle (90 degree)
- Correct speed force
- Ensures stability at needle penetration
- LED light to prevent over-pressing
- Controlled tissue expansion
- Less injection pain
- Less side effects
- Controlled Delivery – Target Deposit
- Aspiration check/avoid direct deposit in vessel
- Needle shield eliminates needle exposure
- Safety lock and Childproof
- Choice of needle extension depth IM/SC Injections
- Keeps needle steady
- Ability to self-inject
- Easy - To - Use

## Injector for convenience and safety

Quality of care for Quality of life

### Specifications & Features

90° Reusable Auto-Injector  
Intramuscular/Subcutaneous  
Dimensions 146 x 39 mm  
Syringes range 0.5 – 5ml  
Needle 30G and larger  
Needle range 8 – 40mm



For more information visit [www.unionmedico.com](http://www.unionmedico.com)

Copyrights all rights reserved by Union Medico - COBIS Science Park - Ole Maaløes Vej 3, 2200 Copenhagen - Denmark